

Nanjing Perlove Medical Equipment Co., Ltd

Add: No.97, Wangxi Road, Jiangning District, Nanjing 211122
Tel: +86-25-87187780 +86-25-68571666-9350 9353
Fax: +86-25-87187780
Http://www.perlove.com.cn
E-mail: overseas@perlove.com.cn




**PERLOVE
MEDICAL**



**OUR UPCOMING
EXHIBITIONS:**

You can see more details on our website
Http://www.perlove.com.cn




PERLOVE MEDICAL

For Overseas Market: Nanjing Perlove Medical Equipment Co.,Ltd.
For Domestic Market: Nanjing Perlove Medical Equipment Co.,Ltd.
Beijing Per-Builders Science and Technology Co.,Ltd.
Zhuhai Perlead Medical Equipment Co.,Ltd.


PLX5100
Mobile Digital Radiography System

Lightweight Design Compact and Flexible

- The rocker arm can move up and down, the combined tube can rotate at multiple angles, with the wireless panel it is more convenient to operate.
- The width of the body is only 62cm, which can easily pass through narrow passages and can be used for clinical exposure in the ward.



- PLX5100 uses lightweight materials, the weight of the whole machine is less than 150KG, innovative mechanical design, light and flexible to move; dislocation brake system, safe and reliable braking.



- Equipped with large size, low energy consumption, portable mobile wireless flat panel detector.
- The workstation adopts DR special image acquisition and processing software, supports DICOM service, convenient to operate.

Mobile DR

Mobile Digital Radiography System



High definition image Easy to clinical diagnosis

- Equiped with high frequency inverter technology to ensure high X-ray quality, low radiation dose and sharp image.

