

DP-50 Vet

Digital Ultrasonic Diagnostic Imaging System

The Premium System for Your Choice



DP-50 Vet

With comprehensive configurations and an integrated design, NEW DP-50 Vet is the outcome of Mindray's continuous and determined efforts towards making veterinary healthcare more efficient, effective and accessible for all. The DP-50 Vet, a premium Black & White ultrasound system, offers more than your expectation. With a smart new shape and full screen, enhanced mobility and more convenient operations, this reliable and affordable system will enhance general practices and is well-suited for multiple veterinary clinical needs.

Enhanced Veterinary Diagnostic Confidence

- **Power Doppler:** Identify blood flow under various veterinary clinical settings
- **PW Doppler and Auto Trace:** Reveal details of blood flow for more comprehensive diagnosis
- **Veterinary Applications:** Measurement, Comment, Report and Body mark designed for canines, felines, equines, bovines, ovines and more



Mobile Design for Veterinary Daily Use

- Sleek, streamlined, and compact shape
- 60 degree tilting angle adjustable monitor
- 2 universal transducer connectors
- 500GB hard disk for large patient data storage
- 1.5h scanning with rechargeable battery
- Slim trolley for easy mobility
- Convenient hand-carried bag for easy transportation

Dedicated Veterinary Workflow

- **iTouch™:** Quick image optimization in one touch
- **iZoom™:** Instant full screen image expansion with one button
- **iStorage™:** Directly transfer images and reports to PC via network cable
- **User-defined keys:** programmable personalized operations tailored for veterinarian



Canine Gallbladder



Canine Bladder

www.mindray.com

P/N:ENG-DP-50Vet-210285x2P-20190918

©2019 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.

 **sonocare**
ARGENTINA

mindray
healthcare within reach